

# STIEFEL BIOPSY PUNCHES diameter 3 mm



Product Code	22651
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

EAN/UPC 5011091105036

## Description

### Stiefel biopsy punches

Disposable sterile instrument to perform an accurate skin biopsy.  
Stainless steel seamless sharp edge.  
Ribbed handle for more precise grip.

Size engraved on each punch for easy identification.

- Stiefel code: STIEF-BP3
- Biopsy punches: Ø 3 mm
- STERILE