

TRUBOND NON ABSORB. SUTURE gauge 2/0 circle 1/2
needle 17 mm - 90 cm - green



Product Code	22770
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2380199
CND	H010201010202
EAN/UPC	8023279227703

Description

TRUBOND - non absorbable surgical suture, silicone coated and braided polyester for permanent support and strength

Equivalent Ethicon brand: ETHIBOND/MERSILENE

- Gauge: 2/0
- Metric: 3
- Colour: Green
- Circle: 1/2
- Needle profile: Taper cutting double arm
- Needle length: 17 mm
- Suture length: 90 cm

Technical details:

- braided polyester, coated with silicone for a smooth knot tie down
- high tensile strength for reliability
- low tissue reactivity - undergoes no significant change in body
- excellent handling and long lasting
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- used for permanent and removable fixtures
- prosthetic implantations
- "Face lifts"
- cardiovascular surgery
- ideal for valve repair/ replacement

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.