

# TRUBOND NON ABSORB. SUTURE gauge 2 circle 1/2 needle 45 mm - 100 cm - green



Product Code	22777
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2380199
CND	H010201010202
EAN/UPC	802327922772

## Description

TRUBOND - non absorbable surgical suture, silicone coated and braided polyester for permanent support and strength

Equivalent Ethicon brand: ETHIBOND/MERSILENE

- Gauge: 2
- Metric: 5
- Colour: Green
- Circle: 1/2
- Needle profile: Taper cutting heavy
- Needle length: 45 mm
- Suture length: 100 cm

### Technical details:

- braided polyester, coated with silicone for a smooth knot tie down
- high tensile strength for reliability
- low tissue reactivity - undergoes no significant change in body
- excellent handling and long lasting
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- used for permanent and removable fixtures
- prosthetic implantations
- "Face lifts"
- cardiovascular surgery
- ideal for valve repair/ replacement

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.