

# TRUSTEEL NON ABSORB. SUTURE gauge 1 circle 1/2 needle 40 mm - 45 cm - natural



Product Code	22790
Unit of sale	box of 48 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	2381046
CND	H0102010199
EAN/UPC	8023279227901

## Description

TRUSTEEL - non absorbable monofilament steel surgical suture

Equivalent Ethicon brand: ETHISTEEL

- Gauge: 1
- Metric: 4
- Colour: Natural
- Circle: 1/2
- Needle profile: Round body blunt point
- Needle length: 40 mm
- Suture length: 45 cm

Technical details:

- monofilament stainless steel LVM 316 grade wire
- optimal compatibility with body tissue
- low or no tissue reactivity
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- orthopedic surgeries
- sternum closure

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.