

# TRULON NON ABSORB. SUTURE gauge 2/0 circle 3/8 needle 26 mm - 45 cm - blue



Product Code	22800
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781306
CND	H010201010102
EAN/UPC	8023279228007
GMDN	16584

## Description

TRULON - polyamide nylon non absorbable monofilament suture  
Equivalent Ethicon brand: ETHILON / NUROLON

- Gauge: 2/0
- Metric: 3.0
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 26 mm
- Suture length: 45 cm

### Technical details:

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.