

TRULON NON ABSORB. SUTURE gauge 6/0 circle 3/8
needle 12 mm - 45 cm - blue



Product Code	22820
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	1781308
CND	H010201010102
EAN/UPC	8023279228205
GMDN	16584

Description

TRULON - polyamide nylon non absorbable monofilament suture
Equivalent Ethicon brand: ETHILON / NUROLON

- Gauge: 6/0
- Metric: 0.7
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 12 mm
- Suture length: 45 cm

Technical details:

- monofilament polyamide
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

Applications:

- for cuticular and subcuticular closure
- orthopedic surgery
- ophthalmic
- general closure
- micro surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.