

# TRUSILK NON ABSORB. SUTURE gauge 3/0 circle 1/2 needle 25 mm - 75 cm - black



Product Code	22842
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781310
CND	H010202020302
EAN/UPC	8023279228427
GMDN	13910

## Description

TRUSILK - braided black silk non absorbable multifilament suture  
Equivalent Ethicon brand: MERSILK / PERMA-HAND SILK

- Gauge: 3/0
- Metric: 2.0
- Colour: Black
- Circle: 1/2
- Needle profile: taper cutting
- Needle length: 25 mm
- Suture length: 75 cm

### Technical details:

- made from 100% protein fiber spun
- coating of wax/silicone
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- general surgery
- ophthalmic surgery
- cardiovascular surgery
- gastrointestinal surgery
- orthopedic surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.