

# TRUSILK NON ABSORB. SUTURE gauge 4/0 circle 3/8 needle 19 mm - 45 cm - black



Product Code	22849
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781312
CND	H010202020302
EAN/UPC	8023279228496
GMDN	13910

## Description

TRUSILK - braided black silk non absorbable multifilament suture  
Equivalent Ethicon brand: MERSILK / PERMA-HAND SILK

- Gauge: 4/0
- Metric: 1.5
- Colour: Black
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 19 mm
- Suture length: 45 cm

### Technical details:

- made from 100% protein fiber spun
- coating of wax/silicone
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- general surgery
- ophthalmic surgery
- cardiovascular surgery
- gastrointestinal surgery
- orthopedic surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.