

# TRULENE NON ABSORB. SUTURE gauge 3/0 circle straight needle 60 mm - 70 cm - blue



Product Code	22875
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No
RDM (NSIS)	1781302
CND	H010201010403
EAN/UPC	8023279228755
GMDN	13909

## Description

TRULENE - monofilament non absorbable suture - polypropylene  
Equivalent Ethicon brand: PROLENE

- Gauge: 3/0
- Metric: 2.0
- Colour: Blue
- Circle: Straight
- Needle profile: conventional cutting
- Needle length: 60 mm
- Suture length: 70 cm

### Technical details:

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.