

# TRULENE NON ABSORB. SUTURE gauge 5/0 circle 3/8 needle 16 mm - 70 cm - blue



Product Code	22881
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781304
CND	H010201010403
EAN/UPC	8023279228816
GMDN	13909

## Description

TRULENE - monofilament non absorbable suture - polypropylene  
Equivalent Ethicon brand: PROLENE

- Gauge: 5/0
- Metric: 1.0
- Colour: Blue
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 16 mm
- Suture length: 70 cm

### Technical details:

- monofilament polypropylene
- non-adhering and biologically inert
- sterilized by ethylene oxide
- shelf life of 5 years

### Applications:

- hernia repair
- permanent fixtures
- general surgery
- micro surgery
- orthopedic surgery
- cardiovascular surgery
- plastic surgery
- can also be used in ophthalmic and neurological procedures

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.