

# TRUBARB ABSORB. SUTURE gauge 3/0 circle 1/2 needle 26mm - 45cm - violet



Product Code	22895
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2380161
CND	H0101010101
EAN/UPC	8023279228953

## Description

TRUBARB - absorbable surgical needed-suture (synthetic) knotless tissue-closure device

Equivalent Ethicon brand: STRATAFIX

- Gauge: 3/0
- Metric: 2
- Colour: Violet
- Circle: 1/2
- Needle profile: Taper point
- Needle length: 26 mm
- Suture length: 45 cm

### Technical details

- monofilament polydioxanone
- secure and effective triangular end-stopper
- single angle cut uni-directional helical barb design
- world class 300 series QNPL needle
- sterilized by ethylene oxide
- shelf life of 3 years

### Applications

- gynaecological and urology procedures
- general surgery procedures
- oncology procedures
- plastic and cosmetic surgeries

Multilingual manual and box: GB, FR, IT, ES, PT, DE, GR.