

MONOGLYDE ABSORB. SUTURE gauge 3/0 circle 3/8 needle 25mm - 70cm - undyed



Product Code	22903
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No
RDM (NSIS)	1781296
CND	H010101010302
EAN/UPC	8023279229035
GMDN	45814

Description

MONOGLYDE - monofilament absorbable suture poliglecaprone 25
Equivalent Ethicon brand: MONOCRYL

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 3/8
- Needle profile: slim blade cutting
- Needle length: 25 mm
- Suture length: 70 cm

Technical details:

- monofilament
- poliglecaprone 25
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 3 years
- wound support: 21-28 days
- mass absorption: 90-120 days

Applications:

- used for subcuticular repair and general soft-tissue approximation
- bowel anastomoses
- closure of the peritoneum
- closure of uterus and vaginal cuff
- not recommended for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery
- since this is an absorbable suture, it should not be used for extended approximation of tissue (fascia)

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.