

# TRUSYNTH FAST ABSORB. SUTURE gauge 3/0 circle 1/2 needle 17mm - 45cm - undyed



Product Code	22911
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781317
CND	H010101010202
EAN/UPC	8023279229110
GMDN	17471

## Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast  
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 1/2
- Needle profile: round body
- Needle length: 17 mm
- Suture length: 45 cm

### Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

### Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.