

TRUSYNTH FAST ABSORB. SUTURE gauge 3/0 circle 1/2 needle 17mm - 70cm - undyed



| | |
|---------------|----------------|
| Product Code | 22913 |
| Unit of sale | box of 12 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | III |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 1781317 |
| CND | H010101010202 |
| EAN/UPC | 8023279229134 |
| GMDN | 17471 |

Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 1/2
- Needle profile: round body
- Needle length: 17 mm
- Suture length: 70 cm

Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.