

# TRUSYNTH FAST ABSORB. SUTURE gauge 3/0 circle 3/8 needle 24mm - 70cm - undyed



|               |                |
|---------------|----------------|
| Product Code  | 22916          |
| Unit of sale  | box of 12 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | III            |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 1781318       |
| CND        | H010101010202 |
| EAN/UPC    | 8023279229165 |
| GMDN       | 17471         |

## Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast  
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 3/0
- Metric: 2.0
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 24 mm
- Suture length: 70 cm

### Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

### Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.