

# TRUSYNTH FAST ABSORB. SUTURE gauge 5/0 circle 3/8 needle 16mm - 70cm - undyed



Product Code	22919
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781318
CND	H010101010202
EAN/UPC	8023279229196
GMDN	17471

## Description

TRUSYNTH FAST - braided & coated polyglactin 910 fast  
Equivalent Ethicon brand: VICRYL RAPID

- Gauge: 5/0
- Metric: 1.0
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 16 mm
- Suture length: 70 cm

### Technical details:

- synthetic absorbable surgical suture
- braided and coated polyglactin 910
- sterilized by ethylene oxide
- shelf life of 5 years
- absorption period: 42 days

### Applications:

- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.