

PDSYNTH ABSORB. SUTURE gauge 0 circle 1/2 needle 40mm - 70cm - violet



Product Code	22942
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781297
CND	H0101010103
EAN/UPC	8023279229424
GMDN	16584

Description

PDSYNTH - monofilament absorbable suture polydioxanone
Equivalent Ethicon brand: PDS or PDS II

- Gauge: 0
- Metric: 3,5
- Colour: Violet
- Circle: 1/2
- Needle profile: round body
- Needle length: 40 mm
- Suture length: 70 cm

Technical details:

- monofilament
- polydioxanone
- synthetic absorbable suture
- sterilized by ethylene oxide
- shelf life of 5 years
- wound Support: 35-42 days
- mass absorption: 180-210 days

Applications:

- for rectus sheath closure
- fascia closure and gastrointestinal tract
- pediatric cardiovascular surgery
- geriatric surgeries
- orthopedic surgery
- general surgery

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.