

TRUSYNTH ABSORB. SUTURE gauge 2/0 circle 1/2 needle 30mm - 90cm - violet



Product Code	22953
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781315
CND	H010101010202
EAN/UPC	8023279229530
GMDN	17471

Description

TRUSYNTH - braided & coated polyglactin 910
Equivalent Ethicon brand: VICRYL

- Gauge: 2/0
- Metric: 3.0
- Colour: Violet
- Circle: 1/2
- Needle profile: round body
- Needle length: 30 mm
- Suture length: 90 cm

Technical details:

- braided and coated polyglactin 910
- coating of calcium stearate and polyglactin 370
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 56-70 days

Applications:

- ophthalmic procedures
- mucosal closures
- skin closure
- dental surgery
- orthopedic surgery
- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.