

# TRUSYNTH ABSORB. SUTURE gauge 4/0 circle 1/2 needle 19mm - 75cm - violet



Product Code	22955
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781315
CND	H010101010202
EAN/UPC	8023279229554
GMDN	17471

## Description

TRUSYNTH - braided & coated polyglactin 910  
Equivalent Ethicon brand: VICRYL

- Gauge: 4/0
- Metric: 1.5
- Colour: Violet
- Circle: 1/2
- Needle profile: taper cutting
- Needle length: 19 mm
- Suture length: 75 cm

### Technical details:

- braided and coated polyglactin 910
- coating of calcium stearate and polyglactin 370
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 28-30 days
- mass absorption: 56-70 days

### Applications:

- ophthalmic procedures
- mucosal closures
- skin closure
- dental surgery
- orthopedic surgery
- general soft tissue approximation

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.