

ETHICON ETHILON MONOFILAMENT SUTURES - gauge 3/0 needle 19 mm PS-2



Product Code	22296
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561746
CND	H010201010101
EAN/UPC	8023279222968
GTIN	30705031001873
GMDN	38000

Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 1669H
- Gauge: 3/0
- Metric: 2.0
- Colour: black
- Circle: 3/8
- Needle length: 19 mm
- Needle: PS-2
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Ethilon is indicated for use in general soft tissue approximation and/or ligation, including cardiovascular, ophthalmic and neurosurgical procedures.