

# TRUGLYDE ABSORB. SUTURE gauge 5/0 circle 1/2 needle 17mm - 75cm - undyed



Product Code	22972
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1781301
CND	H010101020103
EAN/UPC	8023279229721
GMDN	13908

## Description

TRUGLYDE - braided & coated polyglycolic acid

- Gauge: 5/0
- Metric: 1.0
- Colour: Undyed
- Circle: 1/2
- Needle profile: round body
- Needle length: 17 mm
- Suture length: 75 cm

Technical details:

- braided and coated polyglycolic acid
- coating of calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound Support : 28-30 days
- mass Absorption: 70-90 days

Applications:

- for soft tissue approximation during general, gynecological and orthopedic surgeries

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.