

TRUGLYDE FAST ABSORB. SUTURE gauge 1 circle 1/2 needle 40mm - 90cm - undyed



Product Code	22980
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1920380
CND	H010101020103
EAN/UPC	8023279229806
GMDN	13908

Description

TRUGLYDE FAST - Braided and coated Polyglycolic Acid fast Undyed

- Gauge: 1
- Metric: 4.0
- Colour: Undyed
- Circle: 1/2
- Needle profile: round body
- Needle length: 40 mm
- Suture length: 90 cm

Technical details:

- braided and coated
- coated with calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 7-10 days
- mass absorption: 40-42 days

Applications:

- for episiotomy
- for circumcision
- closure of skin

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.