

# TRUGLYDE FAST ABSORB. SUTURE gauge 4/0 circle 3/8 needle 20mm - 75cm - undyed



Product Code	22988
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1920389
CND	H010101020103
EAN/UPC	8023279229882
GMDN	13908

## Description

TRUGLYDE FAST - Braided and coated Polyglycolic Acid fast Undyed

- Gauge: 4/0
- Metric: 1.5
- Colour: Undyed
- Circle: 3/8
- Needle profile: reverse cutting
- Needle length: 20 mm
- Suture length: 75 cm

Technical details:

- braided and coated
- coated with calcium stearate and polycaprolactone
- sterilized by ethylene oxide
- shelf life of 5 years
- wound support: 7-10 days
- mass absorption: 40-42 days

Applications:

- for episiotomy
- for circumcision
- closure of skin

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.