

# ETHICON ETHILON MONOFILAMENT SUTURES - gauge 6/0 needle 16 mm



|               |                |
|---------------|----------------|
| Product Code  | 22300          |
| Unit of sale  | box of 36 pcs  |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | III            |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2561875       |
| CND        | H010201010101 |
| EAN/UPC    | 8023279223002 |
| GMDN       | 38000         |

## Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 660H
- Gauge: 6/0
- Metric: 0.7
- Colour: black
- Circle: 3/8
- Needle length: 16 mm
- Needle: FS-3
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.