

CLINICEL FIBRIL TYPE REGENERATED CELLULOSE

2.5x5.1 cm



Product Code	23007
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1920418
CND	M040501
EAN/UPC	8023279230079
GMDN	38771

Description

CLINICEL OXIDIZED REGENERATED CELLULOSE
Equivalent Ethicon brand: SURGICEL

- Clinicel fibril type
- 2.5 x 5.1 cm
- 1 x 2 inches

Features:

- absorbable hemostat
- made from oxidized regenerated cellulose (plant based product)
- malleable and flexible

Applications:

- Neurosurgery: operative intervention into the cerebrum, ventricular system, spinal surgery
- Cardiac Surgery: valve replacement surgery, sterna medullary bleeding, approaches to the skull base, cerebrovascular surgery
- Gynecological surgeries: uterine wall perforation
- ENT: posterior epistaxis
- General surgeries: iatrogenic perforations

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.

Technical Specifications

Effective Hemostasis: in 3-5 minutes
Sterilization: gamma radiation
Shelf life: 3 years