

# TRULENE NON ABSORBABLE MESH 7.6x15cm - transparent



Product Code	23013
Unit of sale	box of 6 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	1876149
CND	P900202
EAN/UPC	8023279230130
GMDN	60300

## Description

TRULENE (Polypropylene) MESH - Non-Absorbable  
Equivalent Ethicon brand: PROLENE MESH

- Size: 7.6 x 15 cm

### Features:

- non-adhering to tissue and biologically inert
- gives strong support for indefinite periods
- excellent tensile strength and pliability
- easy to draw through tissues
- inert material
- does not harbor micro-organisms
- least thrombogenic

### Applications:

- hernia repair, plastic, general and gastrointestinal surgery

### Technical details:

- monofilament polypropylene non-absorbable surgical sutures
- the product does not require any coating is sold untreated
- sterilized by ethylene oxide

Multilanguage manual and box: GB, FR, IT, ES, PT, DE, GR.

Shelf life of 5 years.