

ETHICON ETHILON MONOFILAMENT SUTURES - gauge 4/0 needle 19 mm



Product Code	22302
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561753
CND	H010201010101
EAN/UPC	8023279223026
GTIN	30705031239542
GMDN	38000

Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 662SLH
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2S
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.