

SUTURE TRAINING PAD WITH WOUNDS WITH MESH



Product Code 23020

Unit of sale 1 pc

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279230208

Description

Silicone suture training pad with 14 wounds.

3 layers: skin, subcutis, muscle (extra layer mesh under the skin).

Soft true to life skin texture.

Multilingual label and manual: GB, FR, IT, ES, PT.

Size: 17.2 x 12 x 1.2 cm