



Product Code 23022

Unit of sale 1 kit

Minimum order 1

Type No medical device

Class

UK-REP Not required

CH-REP Not required

EAN/UPC 8023279230222

Description

SUTURE TRAINING KIT

Suture training kit includes: silicone pad (23020), 6 tools and 4 blades made of stainless steel, 8 sutures and a black PVC leather case.

Intended for teaching and training only, not for human or vet use.

Multilingual label and manual: GB, FR, IT, ES, PT, DE, SE, DK, GR, RO.

Kit size 19.5x22.5x3.5 cm

Standard accessories

1 Case
1 Suture pad

6 Instruments:

1 Needle holder - 12.5 cm
1 Adson forcep - 12.5 cm
1 Iris scissor - 12.5 cm
1 Scalpel handle n.3 - 12.5 cm
1 Mosquito forcep - 12.5 cm
1 Dressing forcep - 12.5 cm

4 Blades:

2 Scalpel blades n.10
2 Scalpel blades n.11

8 Sutures:

2 Nylon sutures 3/0
2 Nylon sutures 4/0
2 Silk sutures 3/0
2 Silk sutures 4/0