

# ETHICON PERMA-HAND SILK SUTURES - gauge 4/0 needle 17 mm J-1- braided



Product Code	22310
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2471021
CND	H010202020301
GTIN	30705031088508
GMDN	13910

## Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 18501G
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 1/2
- Needle length: 17 mm
- Needle: J-1
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.