

# ETHICON PERMA-HAND SILK SUTURES - gauge 3/0 needle 22 mm - braided



Product Code	22311
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2471081
CND	H010202020301
PARAF	964887158
GTIN	30705031007066
GMDN	13910

## Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 632G
- Gauge: 3/0
- Metric: 2.0
- Colour: black
- Circle: 1/2
- Needle length: 22 mm
- Needle: X-1
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.