

ETHICON PERMA-HAND SILK SUTURES - gauge 2/0 needle 22 mm - braided



Product Code	22312
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2471086
CND	H010202020301
PARAF	964887145
GTIN	30705031009435
GMDN	13910

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 737G
- Gauge: 2/0
- Metric: 3.0
- Colour: black
- Circle: 1/2
- Needle length: 22 mm
- Needle: X-1
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.