

ETHICON PROLENE BLUE MONOFILAMENT SUTURES

- gauge 6/0 needle 16 mm



Product Code	22322
Unit of sale	box of 24 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
EAN/UPC	8023279223224
GMDN	13909

Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: W8873T
- Gauge: 6/0
- Metric: 0.7
- Colour: blue
- Circle: 3/8
- Needle lenght: 16 mm
- Needle: PC-3 Prime
- Suture lenght: 45 cm
- Minimum order: box of 24 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.