



Product Code	23495
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
RDM (NSIS)	1842370
CND	Z1203020501
EAN/UPC	6930251800859
UMDNS	16-157

Description

iHEALTH NEO BP5S ARM B.P. MONITOR - with display

- medically validated blood pressure monitors
- easy to read wide Backlit screen
- measurement of blood pressure (systolic, diastolic) and of heart rate (number of beats per minute)
- detection of arrhythmia
- multiple users and unlimited memories
- view historical readings and WHO classification

Multilingual box and manual: GB, FR, IT, ES, DE, GR, PT, NL.

Technical Specifications

Connectivity: Bluetooth 4.0
 Dimensions: 135x55x12 mm
 Screen dimensions: 30.4x27.8 mm
 Cuff circumference: 22-42 cm
 Measuring method: Oscillometric, with automatic inflation and deflation
 Power: DC: 5 V, 1.0 A, 3.7 V, Li-ion battery
 Range:
 Sys: 60-260 mmHg \pm 3 mmHg
 Dia: 40-199 mmHg \pm 3 mmHg
 Pulse: 40-180 bts/min