



Product Code	23499
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
RDM (NSIS)	1439720
CND	Z1203020501
EAN/UPC	0856362005005
UMDNS	16-157
PARAF	972033625

Description

iHEALTH TRACK CONNECTED ARM B.P.M.

- clinically validated blood pressure (ESH)
- systolic, diastolic and heart rate (number of beats per minute)
- evaluation of measurements as recommended by the WHO
- detection of arrhythmia
- colour changing display
- easy to use and comfortable
- 60 memories
- one-Click: save all your measurements in MyVitals App account

Multilingual box and manual: GB, FR, IT, ES, DE, GR, PT, NL.

Technical Specifications

Connectivity: Bluetooth 4.0
 Size: 98x98x46 mm
 Cuff circumference: 22-42 cm
 Weight: 438 g (cuff + batteries included)
 Measuring method: Oscillometric, with automatic inflation and deflation
 Power: 4x1.5V size AAA batteries
 Battery life: 250 measurements
 Range: Sys: 60-260 mmHg, Dia: 40-199 mmHg
 Pulse rate range: 40-180 bts/min