

iHEALTH VIEW BP7S BLOOD PRESSURE MONITOR - wrist with display



Product Code	23502
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

CH-REP	No
RDM (NSIS)	1528328
CND	Z1203020501
EAN/UPC	0855111003941
UMDNS	16-157

Description

iHEALTH VIEW BP7s WRIST B.P. MONITOR - with display

- medically validated blood pressure monitor
- easy to read wide Backlit screen
- measurement of blood pressure (systolic, diastolic) and of heart rate (number of beats per minute)
- detection of arrhythmia
- best wrist position sensor for optimum measurement
- evaluation of measurements as recommended by the WHO (World Health Organization)

Multilingual box and manual: GB, FR, IT, ES, DE, GR, PT, NL.

Motions Sensor Technology: does not allow to take BP in case BP monitor is not correctly positioned

Technical Specifications

Connectivity: Bluetooth 4.0
Dimensions: 72x74x17.6 mm
Screen dimensions: 44x26 mm
Cuff circumference: 13.5-22 cm
Measuring method: Oscillometric, with automatic inflation and deflation
Power: DC: 5 V, 1.0 A, 1x3.7 V, Li-ion 400 mAh
Range: Sys: 60-260 mmHg
Dia: 40-199 mmHg
Pulse: 40-180 bts/min
Weight: 120 g