

iHEALTH ALIGN BG1 BUTTON SIZE GLUCOSE MONITOR mobile transmission



Product Code	23509
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

CH-REP	No
RDM (NSIS)	1439834
CND	W0201060102
EAN/UPC	0855111003705
PARAF	925944478

Description

iHEALTH ALIGN - BG1 GLUCOMETER

Economic, button size model with advanced technology and good reliability.
A convenient way to measure your blood glucose monitor and transmit by your Phone.

Supplied with 10 lancets, lancing device, 2 lithium batteries CR1620, 10 smartphone protection sleeves and pouch.

iHealth Gluco-Smart App, available in GB, FR, IT, ES, DE, CZ, SK, HR, TR, IR, CN, Arabic.
It allows to manage all your smart iHealth devices and track relations among parameters.

Technical Specifications

Machine size: 52x30x9.5 mm
Measuring method: Amperometric technology using glucose dehydrogenase (FAD-GDH)
Result range: 20 mg/dL~600 mg/dL (1.1 mmol/L~33.3 mmol/L)
Power source: DC 3.0, CR 1620
Blood source: Fresh capillary whole blood
Blood volume: Min. 0.5 microliter