


ETHICON VICRYL PLUS ABSORBABLE SUTURES - gauge 4/0 needle 19 mm - braided



Product Code	22360
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2487103
CND	H010101040201
EAN/UPC	8023279223606
GMDN	45401

Description

VICRYL PLUS BRAIDED - absorbable sterile surgical sutures

Complete absorption in 56-70 days

- Ethicon code: VCP392ZH
- Gauge: 4/0
- Metric: 1.5
- Colour: violet
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 45 cm
- Minimum order: box of 36 pieces

vicryl plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular.