

# ETHICON VICRYL PLUS ABSORBABLE SUTURES - gauge 3/0 needle 17 mm - braided



Product Code	22361
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2487059
CND	H010101040201
GTIN	30705031077441
GMDN	45401

## Description

VICRYL PLUS BRAIDED - absorbable sterile surgical sutures

Complete absorption in 56-70 days

- Ethicon code: VCP305H
- Gauge: 3/0
- Metric: 2.0
- Colour: violet
- Circle: 1/2
- Needle length: 17 mm
- Needle: RB-1 PLus
- Suture length: 70 cm
- Minimum order: box of 36 pieces

vicryl plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular.