

STOCKINETTE LEGGINGS 25x75 cm - medium



Product Code	23659
Unit of sale	box of 60 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2119994
CND	T0299
EAN/UPC	8023279236590
GMDN	35531

Description

STOCKINETTE LEGGING 25x75 cm - medium - sterile - box of 60

For orthopedics and gynaecology

Multilingual instructions and label: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, CE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, Arabic

Made in Italy.