

# B BRAUN INTROCAN SAFETY 3 PUR IV CATHETER 22G 25mm - sterile



Product Code	23664
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
--------	-----

RDM (NSIS)	505109
CND	C0101010202
EAN/UPC	4046964047152

## Description

### BBRAUN I.V. CATHETERS - RADIOPAQUE

INTROCAN SAFETY® - sterile

PUR material softens at body temperature and is proven to have a positive impact on indwell time and phlebitis risk

- INTROCAN SAFETY® PUR
- B. Braun code: 4251128-01
- Gauge: 22 G
- Gauge mm: 0.9x25
- Colour: Blue
- Minimum order: box of 50

Latex, PVC, DEHP free

Multi-Access Septum is designed to control the flow of blood from the catheter hub during needle withdrawal and while exchanging subsequent Luer connections.

Automatic safety device to protect clinicians.

Unique geometry of the Universal Back Cut Bevel aids in accessing difficult veins by providing a highly flexible pathway for easy catheter insertion with less tearing.

Integrated Stabilization Platform is designed to minimize catheter movement to help reduce catheter-related complications.

Double Flashback technology helps to support first stick success through quick visual confirmation that both - needle and catheter - are successfully in the vein.

With multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, PL, HU, RO, CZ, SK, BG, GR, HR, RS, SI, EE, LV, LT, SE,FI, NO, RU, TR.