

MESOTHERAPY LUER NEEDLES 30G 0,30x4 mm - yellow



Product Code	23672
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2648227
CND	A019006
EAN/UPC	8023279236729
PARAF	950435697

Description

MESOTHERAPY LUER NEEDLES - STERILE

- Gauge 30G (Ø 0.30 mm)
- Length: 4 mm
- Colour: yellow
- Box of 100 pcs.
- Sterile

Meso-relle disposable hypodermic needles are multiple sharpened to facilitate penetration. Stainless steel AISI 304 cannula and polypropilene connection for syringe.

Sterilized with ethylene oxide.
Sold in individual blister made of medical paper.

Made in Japan.