



Product Code	23686
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2680434
CND	A019006
EAN/UPC	8023279236866
PARAF	963992767

Description

FILLER LUER NEEDLES - STERILE

- Gauge 27G (Ø 0.40 mm)
- Length: 40 mm
- Colour: grey
- Box of 100 pcs.
- Sterile

Meso-relle disposable hypodermic needles are multiple sharpened to facilitate penetration. Stainless steel AISI 304 cannula and polypropilene connection for syringe.

Sterilized with ethylene oxide.
Sold in individual blister made of medical paper.

Made in Japan.