



Product Code	23700
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2100296
CND	A019006
EAN/UPC	8023279237009

Description

MESO-RELLE MULTI-INJECTOR CONNECTORS

- Type: Straight
- 3 Way-out
- Bout: Luer

Disposable mesotherapy multi-injectors connectors sterilized with ethylene oxide.

Packed in individual, transparent blister, colourless and transparent to enable full examination of the inner chamber. Perfectly resistant to high pressure.