



Product Code	23701
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	2100302
CND	A019006
EAN/UPC	8023279237016

Description

MESO-RELLE MULTI-INJECTOR CONNECTORS

- Type: Straight
- 5 Way-out
- Taper: Luer

Disposable mesotherapy multi-injectors connectors sterilized with ethylene oxide.

Packed in individual, transparent blister, colourless and transparent to enable full examination of the inner chamber. Perfectly resistant to high pressure.