

BD VENFLON IV CATHETERS 20G 32 mm - 391892 - sterile



Product Code	23716
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	35887
CND	C0101010101
PARAF	987474727
GTIN	30382903914525
GMDN	64574

Description

BD VENFLON™ I IV CATHETERS WITH BD INSTAFLASH

- BD code: 391892
- Gauge: 20G
- Gauge mm: 32 mm (1¼")
- Flow rate: 54
- Colour: pink
- Box of 50 pcs.
- Sterile

I.V. Cannula Luer-Lock with a plug attached to the flow control plug. Designed to gain access to peripheral veins of the patient's blood system for rehydration, parenteral nutrition, medication delivery, blood transfusion, and monitoring purposes.

BD Instaflash needle technology allows for immediate visualization of blood along the catheter

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, NL, DK, SE, NO, FI, LT, LV, EE, PL, CZ, SK, HU, SI, HR, RO, BG, GR, UA, TR, RU, KZ, Arabic.