

ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 3/0 needle 19 mm - braided



Product Code	22372
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2436966
CND	H0101010299
EAN/UPC	8023279223729
GTIN	30705031106585
GMDN	17471

Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: V2930G
- Gauge: 3/0
- Metric: 2.0
- Colour: undyed
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 45 cm
- Minimum order: box of 12 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.