

STRAIGHT HUB CONVENTIONAL CATHETER 18G

45mm - sterile



Product Code	23720
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No
RDM (NSIS)	1324199
CND	C0101010102
EAN/UPC	8023279237207
PARAF	987471149
GMDN	40601

Description

I.V. CATHETER WITH STRAIGHT HUB - STERILE

- Conventional no protection
- Gauge: 18G
- Gauge mm: 45 mm (1¾")
- Colour: green
- Box of 50 pcs.
- Sterile

- High flow rate: thin wall kink-resistant polyurethane and FEP radiopaque catheters
- Easy vein entry: ultrasharp V-point introducer needle reduces penetration force by 25 to 30%
- Secure grip and intuitive handling: ribbed housing and bevel-up indicator
- Patient comfort and secure fixation: long hub for generous suturing or taping
- Needle-in-vein aspiration: removable flash plug
- Flast visualization: translucent flash chamber