

SIDEPORT CONVENTIONAL CATHETER 18G 45mm - sterile



Product Code	23730
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1323563
CND	C0101010101
EAN/UPC	8023279237306
PARAF	987471176
GMDN	40601

Description

I.V. CATHETER WITH SIDE PORT - STERILE

- Conventional no protection
- Gauge: 18G
- Gauge: 45 mm (1¾")
- Colour: green

- Hygienic sideport access
- High flow rate
- Easy vein entry
- Secure fixation
- Intermittent therapies
- Stabilization