

SIDEPORT CONVENTIONAL CATHETER 22G 25mm - sterile



Product Code	23732
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	1323567
CND	C0101010101
EAN/UPC	8023279237320
PARAF	987471190
GMDN	40601

Description

I.V. CATHETER WITH SIDE PORT - STERILE

- Conventional no protection
- Gauge: 22G
- Gauge: 25 mm (1")
- Colour: blue

- Hygienic sideport access
- High flow rate
- Easy vein entry
- Secure fixation
- Intermittent therapies
- Stabilization